
**endo**doctor®

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## Content

|                                        |    |
|----------------------------------------|----|
| Name of the product                    | 2  |
| CE Marking                             | 2  |
| Address of the Notified Body           | 2  |
| Address                                | 2  |
| Special information on delivery        | 3  |
| Intended use of the product            | 3  |
| Indications, contraindications         | 4  |
| Warnings                               | 4  |
| Intended users                         | 4  |
| Specifications required for proper use | 5  |
| Pretreatment required before use       | 5  |
| Product Lifetime                       | 5  |
| Notes on preparation                   | 5  |
| Manual preparation                     | 6  |
| Mechanical processing                  | 8  |
| Sterilization                          | 9  |
| Inspection, maintenance and testing    | 9  |
| Checking the mechanics and surfaces    | 10 |
| Common use with other products         | 11 |
| Disposal                               | 13 |
| Serious incidents                      | 13 |

- Fig.1: Individual components of the endoscope  
Fig.2: Cleaning the distal end  
Fig.3: Cleaning the light guide attachment  
Fig.4: Cleaning the Proximal Window  
Fig.5: Connection of a fiber optic cable with thread (compatible to the suppliers Storz, Olympus)  
Fig.6: Connection of a fibre optic cable with snap mechanism (compatible with Wolf). The first sleeve of the fibre optic cable attachment is removed  
Fig.7: To insert the endoscope into the shaft, we make sure that the markings of the shaft and the endoscope match. In these positions the components can be pushed together  
Fig.8: Both components can be firmly connected by turning the locking lever clockwise. Accordingly, the connection is released by turning the locking lever counterclockwise.  
Fig. 9: Placing an endoscopy camera on the standardized eyepiece funnel

Before using your new device for the first time, read this manual carefully and follow the instructions and recommendations. This will help you avoid damage that may result from improper use or incorrect / inadequate preparation.

The instructions for use contain important information necessary for safe, professional and economical operation of these devices. We reserve the right to make future technical changes, so that deviations in the contents or the pictorial representation are possible.

These operating instructions are intended to facilitate the use of endodoctor endoscopes, but they are not instructions for endoscopic procedures.

## Name of the product

Rigid endoscopes / endoscopic accessories

## CE Marking

The CE marking confirms that the product complies with the following directives: Rigid endoscopes from ENDODOCTOR GmbH are manufactured in accordance with DIN EN ISO 13485 and meet the requirements of Regulation (EU) 2017/745 of the European Parliament and of the Council of 5 April 2017.

The CE marking, including the identification number of the notified body, is affixed to the medical device

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## Address

Endodoctor GmbH  
Unterer Damm 15,  
78567 Fridingen, Germany

## Special information on delivery

The endoscope is delivered in a cardboard packaging. The endoscope is additionally padded with foam for protection. Rigid endoscope and accessories must be checked for completeness and possible damage immediately upon receipt. Should the packaging show any damage, please inform Endodoctor GmbH or your local dealer immediately. Transport damages can only be claimed if they are reported immediately after receipt of the goods.

Endoscopes are delivered in non-sterile condition.

Endoscopes are fragile, handle with care.

## Intended use of the product

The endoscope serves to visualize the inside of the patient. It can be viewed both through surgically created accesses and in natural body orifices. Rigid endoscopes can be used in combination with surgical instruments for numerous indications. The manipulation is not done with the endoscope itself, but rather with appropriate accessories such as forceps and scissors, which are inserted through one or more accesses or through a shaft - parallel to the endoscope. The endoscope is not to be used for any other purpose.

The different endoscopic procedures are usually designated by the organ being examined, combined with the suffix "-scopy".

Arthroscopy: Used to visualize joints

Bronchoscopy: Serves to visualize the bronchial tubes

Hysteroscopy: Used to visualize the uterus

Laparoscopy: Used to visualise the abdominal cavity and the organs within it

Otoscopy: Serves to visualize the auditory canal and the eardrum

Sinuscopy: Serves to visualize the maxillary sinuses

Thoracoscopy: Serves to visualize the thoracic cavity and the peritoneum

Urethrosopy: Used to visualize the urethra

Ureterorenoscopy: Used to visualize the ureters and renal pelvis

Cystoscopy: Serves to visualize the bladder

Spine surgery: Serves to display the intervertebral discs

## Indications, contraindications

The use of endoscopes from ENDODOCTOR GmbH is generally indicated for endoscopic applications in diagnosis or minimally invasive surgery. The main goal of endoscopic diagnostics or endoscopic surgery is to protect the tissue.

These instructions for use expressly do not describe or explain any clinical application.

The use of rigid endoscopes is contraindicated if endoscopic procedures are contraindicated. Rigid endoscopes from ENDODOCTOR GmbH are not intended for direct contact with the heart or the Central Nervous System.

## Warnings

The high-power light emitted by the endoscope can cause high temperatures at the end surfaces of the light guide and at the distal end of the endoscope. Burns and tissue damage may result from exposure to high heat at:

- prolonged exposure in narrow lumens and cavities
- immediate contact of the distal end to the tissue

Avoid direct contact with surgical drapes, plastic parts, or other flammable materials when the light source is switched on. There is a risk of fire!

preventive measures:

- Avoid long periods of intensive use with high lighting intensity  
Set the light intensity to the lowest possible range so that the field of view is sufficiently illuminated.  
Never place the fiber optic connector or the distal end of the optic on flammable materials. Never cover the endoscope.
- Avoid contact at the light guide connector and at the distal end of the optics  
Switch off the light source before removing the light cable.
- allow optics and light guide cable to cool down after use

## Intended users

The endoscopes may only be used by trained, medical personnel who have been trained on these products.

## Specifications required for proper use

You have decided to buy an endodoctor endoscope and thus acquired a high-quality product. Nevertheless, this endoscope should be handled with care, as rigid endoscopes are very sensitive to bending stress and impact. Such stresses can lead to serious damage to the endoscope. And are not covered by the warranty. If accessories/components are used, the unrestricted function and purpose of the medical device(s) must be ensured and checked by the user before use.

## Pretreatment required before use


Fig. 1: Individual components of the endoscope

The endoscopes must be dismantled into their individual parts so that all surfaces are accessible to the disinfectant. For this reason, all light guide attachments must be removed before each cleaning, if they are used. In general, three different connection options are available. Thus, there are a maximum of two connections which must be removed for cleaning the endoscope.

## Product Lifetime

Based on a validation of the lifetime of the medical device endoscope, it can be confirmed that it can withstand at least 50 reprocessing cycles including steam sterilization without any deviations from the acceptance criteria. Thus the medical device fulfils all requirements of the specified standards over 50 cycles.

## Notes on preparation

**ATTENTION!** The endoscopes must be cleaned, disinfected and sterilized before each use. This applies in particular to the first use after delivery, as all endoscopes are delivered non-sterile.

**ATTENTION!** It is recommended to reprocess the endoscopes as soon as possible after use.

## Manual preparation

**ATTENTION!** Endoscopes must never be cleaned by ultrasound.

The following procedures for preparation and maintenance are only recommendations and make no claim to completeness. The implementation of the selected preparation procedure is the responsibility of the user with regard to its effectiveness.

| Manual preparation for surgical instruments without joint and without lumen<br>Risk class Semicritical A and Critical A (see RKI/BfArM recommendation) |                                                                                                                                                                                                                                                                                                                                                 |                                                                                                                                                                                                                                                  |
|--------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Preparation steps                                                                                                                                      | Procedure                                                                                                                                                                                                                                                                                                                                       | Notes                                                                                                                                                                                                                                            |
| <b>Warnings</b>                                                                                                                                        | Based on a validation of the lifetime of the medical device endoscope, it can be confirmed that it can withstand at least 50 reprocessing cycles including steam sterilization without any deviations from the acceptance criteria. Defective products must have gone through the entire reprocessing process before being returned for repair. | During manual processing, occupational safety must be observed. e.g. personal protective equipment (PPE)                                                                                                                                         |
| <b>instructions:</b>                                                                                                                                   |                                                                                                                                                                                                                                                                                                                                                 |                                                                                                                                                                                                                                                  |
| <b>Preparation at the place of use</b>                                                                                                                 | Wipe off any coarse dirt from the instruments immediately after use. Do not use fixing agents or hot water (>40°C), as this will fix residues and can influence the cleaning success.                                                                                                                                                           | Fixing agents could be e.g. disinfectants                                                                                                                                                                                                        |
| <b>Transport</b>                                                                                                                                       | Safe storage and closed transport of the instruments to the reprocessing site                                                                                                                                                                                                                                                                   | Industrial safety<br><br>To avoid:<br>• Damage to instruments<br>• Contamination of the environment                                                                                                                                              |
| <b>Preparation before cleaning</b>                                                                                                                     | Not required                                                                                                                                                                                                                                                                                                                                    |                                                                                                                                                                                                                                                  |
| <b>Pre-cleaning:</b>                                                                                                                                   | Rinse under running water (< 40°C)                                                                                                                                                                                                                                                                                                              | For coarse soiling                                                                                                                                                                                                                               |
| <b>Cleaning:</b>                                                                                                                                       | Place in suitable cleaning solution according to the instructions of the manufacturer of the cleaner. All surfaces must be wetted. If necessary, brush in the solution until optical cleanliness is achieved.                                                                                                                                   | Note concentration, temperature, degree of contamination and change intervals of the solution.<br><br>Prevents contamination by splashing<br><br><b>The use of combined cleaning/disinfecting agents does not replace the disinfection step.</b> |
| <b>Ultrasonic bath</b>                                                                                                                                 |  <b>ATTENTION! <u>Do not place endoscopes in the ultrasonic bath</u></b>                                                                                                                                                                                     |                                                                                                                                                                                                                                                  |

|                                         |                                                                                                                                                                                                                                                                                                                                                          |                                                                                                         |
|-----------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------|
| <b>Preparation for disinfection</b>     | <p>Rinse all surfaces until clear water is visible.</p> <p>Drain</p> <p>If contamination remains, clean again</p>                                                                                                                                                                                                                                        | It is important to remove residual pollution and chemicals                                              |
| <b>Disinfection:</b>                    | <p>Place in suitable disinfectant solution according to the instructions of the manufacturer of the disinfectant/medical product. All surfaces must be wetted.</p> <p>Rinsing of all surfaces with drinking water/VE water according to microbiological drinking water quality, taking into account the manufacturer's instructions on rinsing time.</p> | Observe concentration, temperature, effective range, exposure time and change intervals of the solution |
| <b>Drying:</b>                          | Drying with the aid of lint-free disposable fleece.                                                                                                                                                                                                                                                                                                      |                                                                                                         |
| <b>Functional testing, maintenance:</b> | <p>Visual inspection for cleanliness and damage</p> <p>If necessary, repeat the reprocessing process until the instrument is visually clean.</p>                                                                                                                                                                                                         |                                                                                                         |

## Mechanical processing

|                                                                                                                                                                                                                                                                |                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>procedure:</b>                                                                                                                                                                                                                                              | <b>Mechanical processing</b>                                                                                                                                                                                                                                                                                                                                                                                                                               |
| <b>products:</b>                                                                                                                                                                                                                                               | Endoscopes without cavities (working channels).                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>instructions:</b>                                                                                                                                                                                                                                           | Based on a validation of the lifetime of the medical device endoscope, it can be confirmed that it can withstand at least 50 reprocessing cycles including steam sterilization without any deviations from the acceptance criteria. Defective products must undergo the entire reprocessing process before being returned for repair. have gone through.                                                                                                   |
| <b>Reprocessing instructions</b>                                                                                                                                                                                                                               |                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Preparation on site:</b>                                                                                                                                                                                                                                    | Remove coarse dirt from the instruments immediately after use. Do not use fixing agents or hot water (>40°C), as this will fix residues and can affect cleaning success.                                                                                                                                                                                                                                                                                   |
| <b>Transport:</b>                                                                                                                                                                                                                                              | Safe storage in a closed container and transport of the instruments to the reprocessing site to avoid damage of the instruments and contamination to the environment.                                                                                                                                                                                                                                                                                      |
| <b>Prepare for decontamination:</b>                                                                                                                                                                                                                            | The instruments must be disassembled or opened for reprocessing.                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Pre-cleaning:</b>                                                                                                                                                                                                                                           | No special requirements                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| <b>Cleaning:</b>                                                                                                                                                                                                                                               | Place the instruments in an open condition in a tray on the trolley and start the cleaning process<br>1. 1 min. pre-rinse with cold water<br>2. emptying<br>3. 3 min. pre-rinse with cold water<br>4. emptying<br>5. 5 min washing at 55°C with 0,5 % alkaline detergent<br>6. emptying<br>7. 3 min Neutralization with warm tap water (>40°C) and neutralizer<br>8. emptying<br>9. 2 min intermediate rinsing with warm tap water (>40°C)<br>10. emptying |
| <b>Disinfection:</b>                                                                                                                                                                                                                                           | Carry out the mechanical thermal disinfection taking into account the national requirements regarding the A0 value (see EN 15883)                                                                                                                                                                                                                                                                                                                          |
| <b>Drying:</b>                                                                                                                                                                                                                                                 | Drying of the outside of the instruments by the drying cycle of the washer-disinfector. If necessary, additional manual drying can be achieved by using a lint-free cloth. Dry cavities of instruments with sterile compressed air.                                                                                                                                                                                                                        |
| <b>Functional testing, maintenance:</b>                                                                                                                                                                                                                        | Visual inspection for cleanliness; assembly of the instruments, care and function test according to the operating instructions. If necessary, repeat the reprocessing process until the instrument is visually clean.                                                                                                                                                                                                                                      |
| <b>packaging:</b>                                                                                                                                                                                                                                              | Packaging of instruments for sterilization according to ISO 11607 and EN 868                                                                                                                                                                                                                                                                                                                                                                               |
| <b>Sterilization:</b>                                                                                                                                                                                                                                          | Sterilization of the products with fractionated pre-vacuum processes (according to ISO 13060 / ISO17665) in accordance with the respective national requirements.<br>3 pre-vacuum phases with at least 60 millibar pressure<br>Heating to a sterilization temperature of at least 132°C; max 137°C<br>Shortest holding time: 3 min<br>rying time: at least 10 min                                                                                          |
| <b>Storage:</b>                                                                                                                                                                                                                                                | Storage of sterilized instruments in a dry, clean and dust-free environment at moderate temperatures between 5°C and 40°C.                                                                                                                                                                                                                                                                                                                                 |
| <b>Information on the validation of the preparation</b>                                                                                                                                                                                                        | The following test instructions, materials and machines were used in the validation:<br>Detergent: Neodisher Mediclean forte 0.5%,<br>FA Neutraliser: Neodisher Z 0.1%<br>Washer-disinfector: Miele PG 8535, Miele G 7736 CD<br>Steam sterilizer Lautenschläger ZentraCert                                                                                                                                                                                 |
| <b>Additional instructions:</b>                                                                                                                                                                                                                                | If the chemicals and machinery described above are not available, it is the responsibility of the user to validate his process accordingly.                                                                                                                                                                                                                                                                                                                |
| It is the responsibility of the user to ensure that the reprocessing process, including resources, materials and personnel, is suitable to achieve the required results. The state of the art and national laws require that validated processes are followed. |                                                                                                                                                                                                                                                                                                                                                                                                                                                            |

## Sterilization

Autoclavable endoscopes from ENDODOCTOR GmbH are marked with the label "AUTOKLAV"

Sterilisation is carried out using the fractionated pre-vacuum process (EN ISO 17665), taking into account the respective national requirements

- 3 pre-vacuum phases with at least 60 millibar pressure
- Heating to a sterilization temperature of at least 132°C; max. 137°C
  - shortest holding time: 3min
  - drying time: at least 10 min

Warning: Please note that sterilization can only be performed on endoscopes and instruments that are free of all organic materials and cleaning residues. In general, sterilisation is only carried out on clean surfaces.

## Inspection, maintenance and testing

Even if optics are carefully cleaned and disinfected, adhesions can occur on the mechanics and glass surfaces of the optics. In these cases, a cotton swab soaked in 70% isoprophyl alcohol can be used for post-cleaning. Clean with moderate pressure in circular movements.

*Fig. 2: Cleaning the distal end*

*Fig. 3: Cleaning the light guide attachment*

*Fig. 4: Cleaning the Proximal Window*

## Checking the mechanics and surfaces

All surfaces of the endoscope must be undamaged and in particular free of sharp edges.

Pay attention to possible dents, kinks, mechanical/thermal damage.

During a visual inspection of the endoscope, the view should be clear and bright. The edges of the endoscopic image should be sharp and rich in contrast.

Check the endoscope light line by holding the distal end against an external light source (window). The light should be visible at the light guide attachment. The light should correspond to that of the light source in strength and colour.

Never use an endoscopy light source to check the light line. This is too strong due to its radiation and can cause damage to the eye.

## Common use with other products

When connecting to other devices or accessories (e.g. TV adapters, light sources, fibre optic cables, cameras), always follow the instructions for use and safety instructions of the respective manufacturers. ENDODOCTOR endoscopes are manufactured according to the valid ISO 8600 standard. This applies especially to the connections of light guide cables and cameras.

*Fig. 5: Connection of an optical fibre cable with thread  
(compatible to the suppliers \*Storz, \*Olympus)*

*Fig. 6: Connection of a fibre optic cable with snap mechanism (compatible with \*Wolf).  
The first sleeve of the fibre optic cable attachment is removed*

\*Storz, Wolf, Olympus, ACMI are registered trademarks of their respective owners. The names are used to describe the connectivity only.

In many cases endoscopes with corresponding shafts are used. These sheaths are usually used to protect the endoscope, for irrigation and suction or to guide an instrument. For this purpose the shafts are attached to the endoscope by means of a locking mechanism.

*Fig. 7: To insert the endoscope into the shaft, make sure that the markings of the shaft and the endoscope match. In these positions the components can be pushed together*

*Fig. 8: Both components can be firmly connected by turning the locking lever clockwise. Accordingly, the connection is released by turning the locking lever counterclockwise.*

*Fig. 9: Placing an endoscopy camera on the standardized eyepiece funnel*

## Disposal

The implementation of European law in national laws and regulations requires the proper disposal of medical devices. Medical devices must be disposed of separately from household waste. Any return to ENDODOCTOR GmbH may only be made in a cleaned & disinfected condition.

## Serious incidents

Professional operators and users and persons who professionally or commercially or in fulfilment of legal tasks or obligations hand over medical devices to the end user for their own use are obliged according to the provisions of the Medical Device Safety Plan Ordinance (MPSV) to report incidents that have occurred in Germany (under certain conditions also incidents that have occurred in third countries) to the Federal Institute for Drugs and Medical Devices (BfArM) or, according to its responsibilities, to the Paul Ehrlich Institute (PEI).

Incidents are malfunctions, failures or changes in the characteristics or performance or inaccuracies in the labelling or the instructions for use of a medical device which have led, could have led or could lead, directly or indirectly, to death or serious deterioration in the state of health of a patient, user or other person.

The notifications shall be addressed to:

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